

AU - 2004 Consistency Determination Apparatus



Consistency determination apparatus for **quick determination of pulp/stock consistency**. The measured stock is poured on to the grid plate of the equipment and drained under suction.

The filtration is repeated until the filtrate becomes perfectly clear. After drying and weight of the sheet the consistency can be determined.

Extra Optional:

Vacuum pump Rapid Sheet Drier for shortening the drying time, as an optional accessory Couch Plate.

SPECIFICATIONS

Capacity	: 1000 ml
Container	: Top easily removable.
Clamping	: Rapid action & efficient.
Grid Plate	: Fitted with a supporting lattice to facilitate removal of the pulp mat from the wire.

OPTION

Vacuum Pump with electric motor.

Rapid Sheet Drier for shortening the drying time

Digital weighing balance

	Shipping Data	Ordering Data
	Net Weight	Model No.
	Gross weight	Voltage
	Dimension	Frequency

Changes:

UEC products undergo continuous development. The technical data in this catalogue are therefore subject to change